

CSD-SR1D

Floor Standing Surgical Scrub Sink



*Better care
Better comfort!*

CSD-SR1D

FLOOR STANDING SURGICAL SCRUB SINK

DESCRIPTION

CSD-SR1D is a humanized disinfection surgical scrub sink, which provides disinfection and cleaning solutions for a variety of scenarios with it's unique multi-functional faucet, simple and modern shape and high-grade stainless steel.

- FEATURES
- The entire sink is made of high-grade stainless steel, corrosion-resistant, heavy duty.
 - The sensor faucet combines soap, water flow and drying functions.
 - With lighting function.
 - With storage space.

TECHNICAL
PARAMETER

Overall Size	L800*W680*H1800mm
Sink	Stainless steel
Structure	Stainless steel


Sensor faucet

The faucet integrates hand soap, water outlet and drying functions, all of which are operated by sensors.


Integrated light

The top of the mirror is equipped with a light for easy use at night.


Storage cabinet

There is storage space at the bottom for storing items and adding soap.


Slope structure

The slope structure drains water quickly, keeping the sink dry and clean. It reduces the sanitary dead corners at the edge of the sink, making cleaning easier and more thorough.



ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US