

Electric Operating Table

Product introduction

YA-10E is an electric operating table designed for all kinds of general operation. It meets all the requirements of thoracic surgery, abdominal surgery, gynecology and obstetrics, ophthalmology and ENT, urology, colon surgery.

The base, support column, accessory fixing device, and accessories are all made of stainless steel, which is easy to clean and anti-rust.

The foam seamless mattress prevents liquid from penetrating and improves patient comfort.




Electric Operating Table

Main details



Main material

The base material is stainless steel304, which is easier to clean, modern and stable.



Brake system

The table can be moved easily or in a stable position thanks to the special brake system.



Emergency “Stop” button

The table is provide with an emergency stop button on the base, help to stop the operating table running in case of emergency situation.



Kidney Bridge

The kidney bridge can be adjusted up to 120mm through the crank.



Power supply

Power with the fuse, and with a master switch,It needs to be pressed when using .



Technical description

and Table specifications

Table specifications

Lowest position of table top

750 mm



Highest position of table top

1000 mm



Back section up motion

70°



Back section down motion

30°



Reverse trendelenburg

0-30°



Trendelenburg

0-25°



Lateral tilt right

0-20°



Lateral tilt left

0-20°



Leg section outward

0-90°



Head section down

0-90°



Kidney elevator mechanism

120 mm



*The page shows the main parameters, Please contact with us for full parameters or more info.

Technical description

and Table specifications

Technical parameter

	PARAMETER
Overall size	2100*550*700-950mm
Reverse trendelenburg	30°
Trendelenburg	25°
Lateral tilt left/right	20°/20°
Head section up/down	55°/90°
Back section up/down	70°/30°
Leg section down	90°
Leg section outwards	90°
Kidney elevator mechanism	120mm
Slide motion	350mm
Loading capacity	200kg

Technical configuration

	PARAMETER
Bed frame	1piece
Mattress	1set
Anesthesia screen	1piece
Shoulder support	1pair
Arm board	1pair
Body support	1pair
Leg support	1pair
Foot plate	1pair
Hand remote	1piece
Fixing clamp	8pieces

*The page shows the main parameters, Please contact with us for full parameters or more info.

ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution.

We started to manufacture surgical table&lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.



YOU TUBE



MEDIK

**Pls scan the QR code
for more info**

info@medikmedical.com | +86 13962231292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US