

# MK-D100HF

Examination Lamp



*Better care  
Better comfort!*

# MK-D100HF

MK-D100HF Exam Light provides a whiter, brighter light than halogen and consumes less energy for improved facility efficiencies—now you can see true tissue color.

## STANDARD SPECIFICATION:

Illuminance	30000Lux0.5m/5,000Lux1m
Colour temperature	4800±200K
Adjust	Adjust light
Total irradiance	380W/m <sup>2</sup>
Size of light field	100mm
Rated power of bulb	3.3/3W
Power supply voltage	220V/50Hz
bulb	LED bulb

## MAIN FEATURES:

- Sustainable energy efficiency with minimal warming
- One of the advantages of LED inspection lights is their compactness
- Usually used during examination or treatment
- Arm 105 cm articulated arm, arm with on/off switch and integrated electronic transformer.
- Ideal cleaning and disinfection conditions
- Single LED lamp head easy to adjust
- With handle and switch button



DETAILS:



**LED LAMP HEAD**

A lamp head with a switch



**ELECTRONICALLY TRANSFORMER**

Transformer with switch and brightness adjustment knob function.



**GOOSENECK TUBE**

Connected to the lamp head is a gooseneck, which can be bent at will.



**HEIGHT ADJUST**

The height of examination lamp can be adjusted by the knob.



**BASE**

Five-claw base with 5 casters, 2 with brakes.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medikmedical.com | +86 199 6269 8550**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US