

# MC-D10

## Gynecological Table



*Better care  
Better comfort!*

# MC-D10

# Gynecological Table

## DESCRIPTION

MC-D10 is an ergonomic and compact design gynecological table, multi-functions with foot pedal control for hands free operation. Stainless steel frame and moulding mattress platform provides the table with long term lifespan and easy for hygiene procedure. Equipped with comprehensive accessories, the MC-D10 gynecological table is an ideal device for baby delivery, women examination and some kinds of surgery in hospital and clinics.

## FEATURES

- Hand controller is easy to operate.
- Unique profile design which is beautiful and elegant.
- Whole stainless steel structure.
- Equipped with hidden sewage basin.
- The table adopts high density sponge and is once foam molding, with no seam.
- Super mute foot wheel and equipped with brake device.
- Armrest swing out  $\geq 90^\circ$  and can slide horizontally.
- Structure and Chassis Stainless steel.

## STANDARD SPECIFICATION

Model	MC-D10
Bed size	1900*600*650-900mm
HI-LOW	650-900mm
Back section is adjustable	-12° to 60°
Seat section is adjustable(option)	0 to 25° option
Trendelenburg	0 to 15°
Reverse-Trendelenburg	0 to 15°
Auxiliary table swing out	0-90°
Loading capacity	200kgs



**POWER SUPPLY**

power with the fuse and with a master switch and also with battery, it needs to be pressed when using.



**EMERGENCY "STOP" BUTTON**

The table is provide with an emergency stop button on the base, help to stop the operating table running in case of emergency situation.



**FOOT SWITCH**

The doctor needs the table to be adjusted when a certain hand operation is being done, the foot pedal switch can free the hand to adjust the table which is a vast convenience.



**FIXED SYSTEM**

There are 2 knobs at the bottom base of the gynecological table, and turn the knob can fix the gynecological table.



**FOOT SECTION**


The foot section can be rotated to one side and removed for easy delivery and exam.



**ANAESTHETIC SCREEN**



**SHOULDER SUPPORT**



**ARMREST**



**LEG SUPPORT**

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medikmedical.com | +86 199 6269 8550**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US