

NEGATOSCOPE

YA-NS01/YA-NS02/YA-NS03



*Better care
Better comfort!*

NEGATOSCOPE

YA-NS01/YA-NS02/YA-NS03

DESCRIPTION

The YA-NS01/02/03 negatoscope with extreme slim body and LED bulb is for all kinds of medical film viewing, it can be applicated for General Film, Digital Film, Breast Mam-mography Film.With several Visual Sizes to choose from ,adjustable brightness, automatic induction lighting.

Series Model	Single YA-NS01	Double YA-NS02	Triple YA-NS03
External Size	478*506*25 mm	838*506*25 mm	1198*506*25 mm
Visual Size	360*425 mm	720*425 mm	1080*425 mm
Max. Power	30W	60W	90W
Light Source	Super High Brighness SMD LED 144 pcs/Panel; 8.000K CT; LED Lifespan of 100,000hrs		
Electric Parameters	Global Voltage Adaptive Built-in Power;AC90V-240V;50/60Hz		
Brightness Setting	3 Digits PWM Dimming System; Brightness 0-100 Stages Continuous Adjust-able among 0-4500cd/m'		
Clip Device	Silicone Clip Device		
Mount Type	Wall-Mounted or Bracket		
Applicable Films	General Film, Digital Film, Breast Mammography Film		
Sensor Function	The Panel Will Be Light Up When Film Inserted and off When Moved away		
Viewing Uniformity	≥90%		
Working Environment	Room Average Illuminance Should No More Than 100 lux		



PWM CONTROLLER FOR DIMMING

The brightness of the light can be adjusted by the knob, Brightness 0-100 level steplessly adjustable 0-4500cd/m.



LED LIGHT SENSOR LIGHTING DESIGN

By pressing the "Auto" button, the mode can be adjusted to sense light, insert the films, and the light will light up, more convenient for doctors to use.



ACRYLIC DIFFUSE PLATE

4mm Imported Acrylic Diffuse Plate, Frosted material, more beautiful, lightweight and tougher at the same time.



OUTSIDE FRAME

25mm thickness, more lightweight, can be hung on the wall or supported by brackets.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US