

# MK-E400M

Field Operation Light



*Better care  
Better comfort!*

# MK-E400M

# FIELD OPERATION LIGHT

## DESCRIPTION

MK-E400M is a classic model field light, it has durable H type base and special detachable design since all the parts can be loaded into the carrying box. It can be installed or dismantled into several parts easily within few mints. Light lux can be adjusted simply when field application.

## FEATURES

- Cold rolled steel with powder coating treatment
- Special knock-down design
- ABS material head unit cover and aluminum bead base
- Lux and spot diameter can be adjust
- Battery back up supply more than 6 hours

## STANDARD SPECIFICATION

Head unit diameter	400 mm
Illuminance (Lux)	> 80,000 lux
Color Temp (K)	4250±250 k
LED Beads qty	36 pcs
LED Beads service life	5,000 hours
Illuminance depth	900 mm
Spot diameter	120-240 mm
Battery back up life	> 6 hours
Left and Right tilting	90 degree
Up and down tilting	Up 35 degree
	Down 50 degree
Carrying box size	800*600*600mm
Packing weight	G.W 45 kg



**Hand push bar**

Here are 3 bars in the head unit  
Knee board can be adjust from 0 to 40 degree.



**Sterilization bar**

Spot diameter can be adjust via the bar.



**Wheels**

Total 4 pcs 3 inch wheels, 2 pcs with locking.



**Power supply**

Power can be used by Inside battery back up or outside power connection.



**Control panel**

Illuminance can be adjust via button.



**FIELD KNOCK-DOWN PACKING**



## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**[www.mediksurgical.com](http://www.mediksurgical.com) | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.