

YA-01F

Field Operating Table



*Better care
Better comfort!*

YA-01F / Field Operating Table

The YA-01F Field operating table based on knock-down design, it is made of light weight high quality raw material only 51.6kg. All the functions can be operated by mechanical operation easily. It will help doctor to adjust to different position shortly and stable.

Main Features

- The table including the accessories can be packed into a small carrying box size only 800*600*600mm and gross weight only 75kg.
- It can be transferred easily and can be installed for using within few minutes when emergency treatment.
- The table base equip with 4 pcs swivel type castors, table can be moved quickly via the wheels.
- Meanwhile there is extra position fix support to guarantee table stability during surgery treatment.



YA-01F / Field Operating Table

Product Details



Manual wheel operation for tilting



Pneumatic control for leg board



Pneumatic control for back board



Ground tolerance adjustment range 0-200mm



Manual wheel operation for Tren and Reverse Tren position



High strength anti-impact material carrying box.
Size 800*600*600mm, G.W 75kg only for easy movement

Technical parameter

	Parameter
Carrying box size	800*600*600mm
Carrying box weight	800*600*600mm
Table size	1980*480*700mm
Table weight	51.6kg
Safe loading weight	135kg
Table height range	700-900mm (±10mm), not include mattress height
Head board adjustment	Up 30°, Down 90°
Back board adjustment	Up 80°, Down 20°
Leg board adjustment	Up 0°, Down 90°, Outward 90°
Tren/Reverse Tren	22°/ 30°
Tilting	24°/ 24°
Castors	100mm, 2 with locking
Mattress	Foam mattress with water-proof cover, thickness 80mm
Ground tolerance adjustment range	0-20mm
Carrying box size	800*600*600mm
Carrying box weight	75kg (include surgical table weight)

Accessories

	Parameter
Anaesthesia screen holder	1pcs
Hand support	1pcs
Body support	2pcs
Leg support	2pcs
Mayo tray	1pcs
IV stand	1 pcs

www.mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.