

# YA-02F

Field Surgery Table



*Better care  
Better comfort!*

# YA-02F

# Field Surgery Table

## DESCRIPTION

YA-02F is a economic light weight design field operation table. It is made of durable steel frame. Suitable in different field application. It includes all basic function to meet different surgery purpose. It will help doctor save more time but can support quick surgery in short time.

## FEATURES

- Cold rolled steel with powder coating treatment
- All operation via mechanical control
- Detachable design for small and quick loading
- Suitable for field and special location application

## STANDARD SPECIFICATION

Size	1980*480*490mm
Height range	490-800mm
Head board	+30 / -90 degree
Back board	-15 / +85 degree
Tren / Reverse Tren position	24 / 24 degree
Ground depth adjustment range	0-110mm
Loading weight	135kg
Carrying box size	1100*650*800mm
Carrying box Weight	75kg

## ACCESSORIES

- Mattress
- Arm support
- Leg support
- Body support
- Carrying box



**Head board adjustment**

Head board can be adjust from +30 to -90 degree via mechanical knobs.



**Castor**

2 pcs dia 100mm heavy duty wheel.



**Support holder**

Ground depth adjustment range can be adjusted via 2 pcs support holder.



**Carrying box**

High strength carrying box size only 1100\*650\*800mm, G.W 75KG.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**[www.mediksurgical.com](http://www.mediksurgical.com) | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US