

# MK-S20

Mayo table



*Better care  
Better comfort!*

# MK-S20

# Mayo table

## DESCRIPTION

MK-S20 is a mayo table made of high grade 304 stainless steel against corrosion, with its unique forklift structure and effortless height adjusting mechanism, it is highly accepted by doctors and nurses in operation room.

## FEATURES

- High grade 304 stainless steel, anti-corrosion and wear-resistant.
- Removable tray ,easy to clean.
- Two castors with brake ,easy to drive and standing.
- long and heavy forklift base design ,more stable.
- Smooth finishing without sharp edge,more safety.
- Height can be adjusted smoothly by gas spring.

## STANDARD SPECIFICATION

Size	600*405*960~1260mm
Structure	Stainless steel
Castor	Dia75mm
Application	Hospital / Clinic
Forklift base length	530mm



**Height adjustment**

The height can be adjusted smoothly by gas spring , which is suitable both for doctor and operation table with different height.



**Removable tray**

Stainless steel removable tray,easy to clean.



**Castors and the forklift base design**

Two castors and both of them with brake which can be stand still in any state , flexible and light to drive ,the long and heavy forklift base will make the table more stable and can be placed under different operating tables.



**Smooth finishing**

All the welds have been carefully polished, the appearance is smooth and textured, and it will be safer to use.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US