

# MK-S21A

Mayo table



*Better care  
Better comfort!*

# MK-S21A

# Mayo table

## DESCRIPTION

MK-S12A mobile height adjustable stainless steel mayo table, all the material is high quality stainless steel , Widely used in the emergency room, nurse station ,operating room and etc.

## FEATURES

- High grade stainless steel
- Smooth finishing without sharp edge
- With four castors, easy to move

## STANDARD SPECIFICATION

Size	600*400*820~1300mm
Structure	Stainless steel
Castor	Dia75mm
Application	Hospital / Clinic



### HEIGHT ADJUST

The height is adjusted by the knob on both sides , which is more stable.



### MATERIAL

We use the high quality stainless steel, Perfectly cut and bent at every corner.



### REMOVABLE TRAY

Stainless steel removable tray,easy to clean.



### CASTORS

Diagonal lockable castors, flexible and light to drive, smooth steering, wear-resistant, two with brake, can be easily braked in any state.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US