

YA-10EB

Surgical Table



*Better care
Better comfort!*

YA-10EB

SURGICAL TABLE

DESCRIPTION

YA-10EB is an electric surgical table which is suitable for the major universal surgeries like thoracic surgery, abdominal surgery, brain surgery, ophthalmology, otolaryngology, obstetrics and gynecology, urology, orthopedics, etc. Stainless steel material and heavy duty structure makes the table durable in use for long lifespan. The table with horizontal sliding function is suitable for specific surgery which needs whole X-ray translucent compatible for C-arm exam.

FEATURES

- 1.Eight independent motors for movements of HI-LO, back section, Trendelenburg and Reverse-Trendelenburg, Lateral tilting function, Horizontal and longitudinal movement(For C-arm examination), Leg section(separated left or right leg section, also left and right leg section as one united set), auto-contour of operating table.
- 2.The table can be moved vertically, and it can be used with C-arm for radiographic diagnosis or filming.
- 3.The countertop is made of high-strength melamine resin sheet, and the mattress is made of memory foam mattress.
- 4.The leg plate can be removed, manually rotated abduction, folded down, easy to adjust, very convenient for urological surgery.
- 5.The handheld remote control which is easy to operate, safe and reliable.
6. Gas spring for adjusting movements of headrest section.
7. Polyurethane shaped mattress which is easy to remove and clean.
8. AISI304 stainless steel structure and chassis, which is easy to clean.
9. Polyurethane shaped mattress which is easy to remove and clean.
10. Gas spring for headrest movement.
11. Mechanism system for waist board lifting Leg section can be removable, then connected to orthopedic traction shelf directly



Remote Controller

The remote control icon is clear and large, which is convenient for quick identification and operation.



Kidney Bridge

The kidney bridge is adjusted by the crank which can up to 110mm.



Form Mattress

High-quality poly surface for easy clean, removable, impermeable.



Stainless Steel Basement

The basement is made of stainless steel 304, which is stable and with a good outlook.



Castors And Lock System

The surgical table is equipped with casters easy to move the surgical table to the place they need , and it with the lock system.

Technical parameter

	PARAMETER
Overall size	2010*500*700-1000mm
Back section up/down	75°/10°
Head section up/down	50°/90°
Leg section down	90°
Reverse trendelenburg	25°
Trendelenburg	25°
Lateral tilt left/right	20°
Slide motion	350mm
Kidney elevator mechanism	110mm
Load capacity	250kg

Back Section Up Motion
-10° to +75°



Leg Section Motion
-90°



Removable Abduction
180°



Slide Motion
350mm



Kidney Elevator Mechanism
120mm



Reverse Trendelenburg
25°



Trendelenburg
25°



Lateral tilt left
20°



Lateral tilt right
20°



info@mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.