

# YA-A602

Lithotomy Stirrup



*Better care  
Better comfort!*

# YA-A602

# LITHOTOMY STIRRUP

## DESCRIPTION

The stirrup provides great stability and comfort for patients during lithotomy or gynecological surgery, can fully expose the patient's surgical area to the doctor. At the same time, the padded and wrapped cushion mattress can avoid iatrogenic injury during long term surgery and the gas assisted adjusting mechanism can reduce the work intensity of the medical staff and improve the operation efficiency.

## FEATURES

- It can meet various body position requirements for lithotomy surgery.
- The pneumatic operating system is light and flexible, reducing the work intensity of medical staff.
- Boots designed to enhance patient comfort and avoid complications.
- The angle can be adjusted precisely to improve the operation efficiency.

## STANDARD SPECIFICATION

Cantilever Length	110cm
Lift up and down	110°
Swing outward and interior	25°
Telescopic distance of sheath	310mm
Leg Weight Capacity	130Kgs
Specifications of side rails of operating table to be applied on	25*10mm /25*8mm /30*10mm /30*8mm /28.5*9.5mm can be customized
Product Weight	13kgs



**Control Handle**

Lift up and down with the handle, easily move the legs to the desired position.



**Fixed Handle**

Loosen the handle to adjust the angle and distance of the leg pedal, and tighten the handle to stabilize the leg pedal in a suitable position.



**Pneumatic Auxiliary Device**

Pneumatic way, with angle adjustment scale, light and stable to use, can stop at any position.



**Boot Design**

The non-compression boot design is adopted, which can fully protect the popliteal fossa, gastrocnemius muscle, peripheral nerves of the hip joint and soft tissues to avoid complications.



**Foldable Design**

Can suit for any surgical table.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**[www.mediksurgical.com](http://www.mediksurgical.com) | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US