

YA-A603

lithotomy stirrup



*Better care
Better comfort!*

YA-A603

lithotomy stirrup

DESCRIPTION

The YA-A603 Lithotomy stirrup is a very practical surgery accessory to improve the patient's comfort and doctor's convenience during surgeries like gynecology, urology, Proctology etc..

The dynamic position adjusting can be achieved for various distance and angles during surgery through easy pneumatic way. And due to the moulded leg fixation system, patient will not feel uncomfot during a long surgery. And this stirrup can almost work compatibly with any brand surgical table.

STANDARD SPECIFICATION

Cantilever Length	100cm
Lift up and down	90°/ -35°
Swing outward and interior	30° / 15°
Telescopic distance of sheath	330mm
Patient Weight Capacity	150Kgs
Specifications of side rails of operating table to be applied on	25*10mm /25*8mm /30*10mm /30*8mm /28.5*9.5mm
Product Weight	13kgs



BOOT DESIGN

The boot design prevents stretching of the groin and lumbar muscles by stopping the knee from falling out in high lithotomy.



LIFT UP AND DOWN

The floating boot reduces the pressure on the calf muscles when the stirrups lift up and down.



REUSABLE PAD

Completely encapsulates the foot, ankle and calf for comfortable and secure positioning.



GERMAN PNEUMATIC AUXILIARY DEVICE

The pneumatic auxiliary device assist reduces the workload for staff and increases safety for all.



CLAMP

Fix lithotomy stirrup on the operating table.



FIXED HANDLE

Adjust the position of lithotomy stirrup , fix lithotomy stirrup.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

www.mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US