

# YA-A505

## Neurosurgery Skull Clamp



*Better care  
Better comfort!*

# YA-A505

# NEUROSURGERY SKULL CLAMP

## DESCRIPTION

The Neurosurgery skull clamp can firmly fix the head in the best surgical position, making sure the surgical site is exposed to the doctor's vision field as much as possible, which is convenient for the doctor to operate. The skull clamp can be adjusted at multiple angles, suitable for various neuro surgeries. The skull clamp can be compatible or customized with any operating bed providing our connection system.

## FEATURES

- Aluminum-magnesium alloy frame and 304 stainless steel accessories, durable and easy to sterilize.
- Surgical site can be exposed to the doctor's field of vision as much as possible.
- The base rod can be customized for various operating tables, multi-angle adjustment, can adapt to various positions of surgery.
- U-shape head clamp and gel pad can suit for different kinds of surgery.

## STANDARD SPECIFICATION

Width adjustment for base rod	80-600mm
Length of U-shape head clamp	500mm
Front and rear adjustment distance	400mm
left and right adjustment distance	450mm
Lifting distance	800mm
Base rod diameter	13-20mm customized



**Storage Suitcase**

All accessories are stored in the suitcase, which is convenient for the doctor to install and also protects the accessories from damage.



**Material**

Aluminum-magnesium alloy plate, the surface is treated with matte hard oxidation to reduce light reflection, wear-resistant and easy to clean and disinfect. High-precision processing makes the parts fit perfectly and provides firm support for surgery.



**Fixing Nails**

Three-point fixation, equipped with two sets of 304 stainless steel nails for adults and children, standard 80lbs nail jacking handwheel with scale, quick locking device, saving operation time.



**Gel Pad**

Equipped with a polymer gel pad, the texture is soft and the elasticity is close to that of human skin, which can relieve local tissue pressure, reduce nerve compression, avoid iatrogenic injury, and can be used for prone position surgery.



**Universal Ball**

The connection is equipped with a universal ball, made of 304 stainless steel, which can be rotated 360 degrees in the plane and 90 degrees longitudinally, suitable for lateral position surgery, equipped with a barrier-free indentation handle, and the connection can be locked at any angle.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**[www.mediksurgical.com](http://www.mediksurgical.com) | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US