

MK-S06

Dressing And Care Cart



MK-S06 DRESSING AND CARE CART

This MK-S06 stainless steel care cart is designed for hospitals, clinics, and nursing environments. Equipped with two storage shelves, a side push handle, and an integrated waste bucket, it provides a practical solution for transporting medical supplies while improving daily nursing efficiency.

Main features

- Premium 304 stainless steel construction for durability and corrosion resistance.
- Two storage shelves with raised edges to help prevent items from slipping.
- Side push handle for comfortable and efficient transportation.
- Integrated 8L stainless steel waste bucket with foot pedal for hands-free operation.
- Four 100mm silent swivel castors, two with brakes for stability.
- Smooth polished surface with rounded edges for easy cleaning and enhanced safety.
- Customization options available upon request.

Main details



Material

High-quality 304 stainless steel construction with excellent corrosion resistance and easy cleaning.



Double Shelf Design

Two spacious stainless steel shelves provide organized storage for medical instruments and nursing supplies.



Side Push Handle

Ergonomic side handle allows smooth and comfortable maneuvering during daily use.



Hands-free Waste Disposal

Equipped with an 8L stainless steel waste bucket operated by foot pedal to help reduce hand contact and improve hygiene during routine nursing procedures.



Castors

Four silent swivel castors ensure smooth mobility, with two locking castors for secure positioning.

Technical parameter

Model	MK-S06
Overall Size	600×400×800mm (23.6" × 15.7" × 31.5")
Material	304 Stainless Steel
Shelves	2 pcs
Waste Bucket	8L Stainless Steel with Foot Pedal
Castors	Dia.100mm (4")
Application	Hospital / Clinic / Nursing Center



ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOUTUBE



MEDIK

Pls scan the QR code for more info

info@mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

©2026 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 01-EN-2026 EN-US