

MK-CS01

Stainless Steel Cabinet



*Better care
Better comfort!*

MK-CS01

STAINLESS STEEL CABINET

DESCRIPTION

This MK-CS01 secure hospital cupboard with glass doors keeps medical supplies organized and visible, improving workflow and saving time.

- FEATURES
- Heavy-duty, all-welded stainless steel construction for durability and strength

- Space-saving stainless steel cabinetry for maximum storage efficiency

- Transparent glass doors for clear visibility of stored items

- Lockable glass doors for enhanced security

- Four adjustable and removable shelves to accommodate various storage needs

TECHNICAL
PARAMETER

Model	MK-CS01
Overall size (L*W*H)	900*400*1780mm
Structure	Stainless steel
Storage shelf	4pcs



Stainless steel handle

With a protruding stainless steel handle on the door for easy operation.



Keyed door latch

Lockable doors provide security for storing a variety of medical supplies.



Shelves

Four removable, adjustable stainless steel shelves.



Transparent Glass Doors

Full height, stainless steel doors with glass panes.



ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US