

# MK-MS04

Mayo table



*Better care  
Better comfort!*

# MK-MS04

## MAYO TABLE

The stable H-shaped base, square tube structure, thick tray and perfect polished finishing allow MK-MS04 to be flexibly used in various surgical situations. The stable height adjustment allows it to be matched with various operating tables, providing doctors with smooth operation and surgical experience.

### FEATURES

- Total SS construction, used mainly in operation theatre
- H-shape tubular base to ensure stability
- Raised edges on all sides of the tray to prevent falling off of objects
- Removable top tray,for easy disinfection and cleaning
- Four antistatic castors with two brakes
- Height adjusted by knob

ITEM	PARAMETER
Overall Size	L600*W400*H750-1100mm
Working shelf	L600*W400mm*L30mm
Safe Working Load	40kg
Height range	H750-1100mm
Castor	Dia75mm
Packing	Knock down design




**Adjustment knob**

With one adjustment knob ,the height of table can be adjusted and stopped at any height.



**Top tray**

The top tray of the mayo table is removable for easy disinfection and cleaning,the raised edge can prevent falling of objects.



**Basement**

The H shape basement with rectangular tube is durable and stable.



**Castors**

4 units silent castors with diagonal lock, easy to move.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**[www.mediksurgical.com](http://www.mediksurgical.com) | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US