

MK-S08D

Medical Dressing Trolley



*Better care
Better comfort!*

MK-S08D

MEDICAL DRESSING TROLLEY

DESCRIPTION

The MK-S08D Dressing Cart is designed for healthcare people to save time and energy by collects dressing materials in one trolley. It keeps dressing materials sterile and reduces the risk of infection with sliding drawer. Additionally, it collects waste materials in a separate place and prevents environmental pollution via the bucket with lid.

- FEATURES
- The trolley is made of stainless steel, which is reliable, smooth and elegant.
 - With two drawer and guardrails around.
 - User-Friendly Handle: A comfortable grip on the top shelf facilitates smooth movement, enhancing usability and convenience.
 - Features two stainless steel drawers, which are effortless to store items.
 - Minimum loading capacity of 200lbs.


STANDARD SPECIFICATION

Model	MK-S08D
Overall size (L*W*H)	670*480*950mm
Structure	Stainless steel
Storage shelf	2pcs
Castors	Four castors with brakes
Drawer	2pc



Solid Steel Construction

With a solid steel structure, the stainless steel surface is easy to clean and not easy to rust.



Detachable waste bin

Adjustable waste bin, the height can be adjusted to satisfy the needs .



Drawers

The drawer with the mute track which can pull out and push very quiet, Large capacity, can put more instruments.



Castors

Equipped with four 360°wheels, all with foot brakes, this rolling cart allows easy transportation of tools or items.



ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



MEDIK

Pls scan the QR code for more info

info@mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US