

MK-S23

Mayo table



*Better care
Better comfort!*

MK-S23

MAYO TABLE

DESCRIPTION

MK-S23 is a hydraulic surgical mayo table, the hydraulic operation system can let the medical staff to adjust the height via foot only with hand free, it can improve the operation experience while minimizing the contamination risk.

FEATURES

- Total SS construction, used mainly in operation theatre.
- U-shape tubular base to ensure stability.
- Raised edges on all sides of the tray to prevent falling off of objects .
- Four antistatic castors with two brakes.
- Lift system with high precision hydraulic pump, pedal-operated.
- Knock down design.

STANDARD SPECIFICATION

Overall Size	L600*W400*H770-1170mm
Working shelf	L600*W400mm
Max. workload	30kg
Height range	770-1170mm
Castor	Dia75mm
G.W.:	16.7kg
Package	Ctn Size: 71.5*52*12cm(1 set/1package)



Foot Pedal

Height adjusted by hydraulic foot pedal



Chassis

U shape chassis,which is stable during transportation



Smooth Finishing

Perfectly polished and welding ,no sharp edge



Tray

Removable tray with raised edge,prevent items from dropping while for easy cleaning



Castors

4 units silent castors with 2 brakes

ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US