

YA-07EH

Electro Hydraulic Operating Table



*Better care
Better comfort!*

YA-07EH ELECTRO HYDRAULIC OPERATING TABLE

DESCRIPTION

The YA-07EH electro-hydraulic operating table has extra-long horizontal sliding together with interchangeable head&foot section that can be suitable for both X-rays and C-arms. It is equipped with an intelligent automatic control system and has variety of flexible body position options to meet the needs of various general surgical operations. The powerful lifting system and extreme load system also make YA-07EH a very welcome option for the hospitals.

FEATURES

- AISI304 stainless steel structure and chassis, which is easy to clean.
- The leg plates adopt imported gas springs, which are easy to control.
- The head and foot section are interchangeable, which is convenient for C-arm or X-ray perspective without blind spots on the bed surface.
- The operating table has one-button reposition function.
- One-button to operate reflex, flex and chair position, easy to operate.
- Maximum 350mm motorised longitudinal shift for flexible use of the C-arm.
- Double control system, hand remote control and Column control.
- Electric brake design, easy to move or fix the operating table.
- Battery charging indicator, beep alarm for lake of power.



Hand remote control

The hand remote control has LCD display. When you press a button, the corresponding icon will be displayed on the LCD display.



Column control

There are also column controls with clearly iconized buttons on the column.



Castor

Square base with castor, 2 diameter 50mm, 2 diameter 120mm, hand remote control brake.



Interchangeable head and foot sections

Modular design allows interchangeable heads and feet.



Kidney Bridge

The built-in kidney bridge can be raised 120mm via the side crank.



Power and switch

Control power, the motor will not work when turned off.

TECHNICAL PARAMETER

Model	YA-07EH
Overall Size(L*W*H)	2160*550*(690-1090)mm
Height Adjustment	690-1090mm
Head section up/down	45°/90°
Back section up/down	90°/40°
Leg section up/down	30°/90°
Leg section outward	90°
Trendelenburg	30°
Reverse trendelenburg	30°
Lateral tilt left/right	20°/20°
Flex&Reflex	90°/220°
Kidney Bridge Up	110mm
Sliding	350mm
Load capacity	320kg

STANDARD ACCESSORY

Name	Quantity
Headrest	1 pc
Armrest	1 pair
Shoulder support holder	1 pair
Waist support holder	1 pair
Leg support holder	1 pair
Footrest	1 pair
Anesthesia screen frame	1 pc
Hand remote control	1 pc
Column control	1 pc

FUNCTIONS

Lowest height position

690mm



Lowest height position

1000mm



Back section up

90°



Reverse trendelenburg

90°



Trendelenburg

90°



Lateral tilt left

20°



Lateral tilt right

20°



Head section up

45°



Head section down

90°



Slide motion

350mm



Leg section up

30°



Head section up



Kidney elevator mechanism

110mm



Chair position



ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US