

YA-B301

Shoulder Attachment



*Better care
Better comfort!*

YA-B301

SHOULDER ATTACHMENT

DESCRIPTION

YA-B301 Shoulder Attachment is specially designed with detachable boards for shoulder surgery which can expose the needed shoulder area to the surgeons clearly. It provides tremendous convenience and comfort for both patients and surgeons for shoulder surgery.

FEATURES

- The bottom bracket is made of 304 stainless steel, and the rest is aluminum-magnesium alloy. The back panel is made of transparent material with waterproof and anti-static high-elastic PU sponge pad.
- The shoulder can be removable, suitable for a variety of shoulder operations.
- It can be adapted to operating tables of various sizes.
- Maintenance-free gas spring power design.

STANDARD SPECIFICATION

| | |
|----------------------------|-------------|
| Backrest size | L560*W520mm |
| Headrest height adjustment | 0-300mm |
| Applicable side rail size | 3mm |



Storage Trolley

Equipped with an auxiliary cart for easy storage and loading and unloading. The base that inserts into the table side rails can be customized to suit a variety of operating tables.



Shoulder Panel Removable

The left and right back panels can be disassembled, so that the surgical site can be exposed to the doctor's field of vision as much as possible.



Backrest Adjustment

The headrest can adjust the front and rear distance, up and down height, and 360-degree universal angle. The inclination angle of the back is adjusted by the gas spring, and the gas spring adopts imported maintenance-free gas spring.



ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

www.mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US