

# YA-CT04

Mortuary Lifting And Transport Trolley



*Better care  
Better comfort!*

# YA-CT04

# MORTUARY LIFTING AND TRANSPORT TROLLEY

## DESCRIPTION

Whole high grade stainless steel material, double motorized lifting mechanism, the features make the mortuary lifting and transport trolley with superior stability, effortless operation and lasting for long-term use.

## MAIN FEATURES

1. The lifting trolley is composed of chassis, lifting bracket moving tray, etc.
2. The whole body is made of stainless steel to prevent rust and corrosion.
3. It is specially for autopsy table or body freezer.
4. The lifting trolley adopts an X-shaped structure, with strong carrying capacity and large stroke.
5. Electric control height stroke 1000mm and easy to operate.
6. Equipped with 6 casters for easy mobility.


## TECHNICAL PARAMETER

Model	YA-CT04
Size	2140*700*400-1400mm
Height adjustment	400-1400mm
Capacity loading	250kg

## PRODUCT DETAILS

### Structure

The overall structure is scissor X-shaped, with large load-bearing capacity and high lifting stroke.




### Lift switch

Control the overall height of the corpse lift table via push button switch




### Motor

Equipped with 2 motors to control height, height stroke 1000mm




### Pallet holder

Equipped with 4 pallet holders to fix the pallet of the corpse lifting table



### Power control box

Equipped with a stainless steel power control box, which will not affect the overall operation and is safe and stable



## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**[www.mediksurgical.com](http://www.mediksurgical.com) | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US