

YA-SR2B

Surgical Scrub Sink



*Better care
Better comfort!*

YA-SR2B

SURGICAL SCRUB SINK

DESCRIPTION

YA-SR2B is an intelligent scrub sink designed in accordance with ergonomic principles. Made of top-grade 304 stainless steel, it has a spacious use space, making it convenient for medical staff to clean and disinfect their hands efficiently.

- FEATURES
- All stainless steel construction

-All inside corners are rounded for ease in cleaning

- Unique slope splash-proof structure

- Infrared sensor control or knee control

-Humanized design

TECHNICAL
PARAMETER

Material	High-quality stainless steel
Size(L*H*T)	1500*820*700mm
304 stainless steel thickness	1.2mm
Options	Quick-release base
Control method	Infrared sensor control or knee control




Sensor faucet

Induction water, water temperature can be adjusted freely.



Arc-shaped groove

The curved trough design brings medical staff closer to the faucet, improving the comfort and convenience of use.




Integrated LED light mirror

Convenient lighting, high-definition mirror.



Slope structure

The slope structure drains water quickly, keeping the sink dry and clean. It reduces the sanitary dead corners at the edge of the sink, making cleaning easier and more thorough.



ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US