

YA-SR2E

Wall-mounted Double
Station Scrub Sink



*Better care
Better comfort!*

YA-SR2E

WALL-MOUNTED DOUBLE STATION SCRUB SINK

DESCRIPTION

YA-SR2E is a double-station wall-mounted scrub sink with a simple stainless steel appearance. It is equipped with induction and knee-controlled water outlets. The water outlet sensing distance is adjustable and it is equipped with a touch-screen operation panel. The water outlet temperature can be set individually. It has a humanized design and is splash-proof and corrosion-resistant.

- FEATURES
- Morden design and easy to clean.
 - Sensor faucet and knee control design.
 - Stainless steel structure, heavy duty.
 - With touch panel, the water temperature of each tap can be adjusted individually.
 - Optional induction soap box, integrated LED light mirror, mechanical or electronic constant temperature water mixing regulating valve.

TECHNICAL
PARAMETER

Overall Size	L1500*H850*D650mm
Sink	Stainless steel
Structure	Stainless steel



Sensor faucet

Induction water, water temperature can be adjusted freely.



Touch panel

The touch panel can adjust the water temperature separately, humanized design.



Knee control design

Knee control panel is more hygienic and convenient.



Integrated LED light mirror

Convenient lighting, high-definition mirror.



Slope structure

The slope structure drains water quickly, keeping the sink dry and clean. It reduces the sanitary dead corners at the edge of the sink, making cleaning easier and more thorough.

ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US