

# YA-SR2F

Double Station Scrub Sink



*Better care  
Better comfort!*

YA-SR2F

DOUBLE STATION SCRUB SINK

DESCRIPTION

YA-SR2F is a two-station scrub sink with a stainless steel body that is corrosion-resistant and antibacterial. It is equipped with an induction faucet and knee-controlled design, which is more hygienic and easy to operate.

- FEATURES
- Morden appearance ,easy to clean
  - Sensor faucet and knee control design
  - Stainless steel structure, heavy duty
  - Optional induction soap box, integrated LED light mirror, mechanical or electronic constant temperature water mixing regulating valve.

TECHNICAL  
PARAMETER

Overall Size	L1500*H850*D650mm
Sink	Stainless steel
Structure	Stainless steel




**Sensor faucet**

Induction water, water temperature can be adjusted freely.



**Knee control design**

Knee control panel is more hygienic and convenient.



**Integrated LED light mirror**

Convenient lighting, high-definition mirror.



**Slope structure**

The slope structure drains water quickly, keeping the sink dry and clean. It reduces the sanitary dead corners at the edge of the sink, making cleaning easier and more thorough.



## ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@mediksurgical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US