YA-02E Plus

Electric Ophthalmic Operation Stretcher



MEDIK: SURGICAL ROOM & NURSING CARE SOLUTIONS
+86 139 6223 1292 | www.mediksurgical.com

YA-02E PLUSELECTRIC OPHTHALMIC OPERATION STRETCHER

The YA-02E Plus electric ophthalmic operation stretcher is a multifunctional, mobile surgical table designed specifically for ophthalmic procedures. It allows for smooth adjustments to height, head rest, backrest, and leg adjustment, accommodating the various patient positioning requirements of ophthalmic and ENT surgeries. the central wheel mechanism facilitates movement within the operating room. With its compact design and easy-to-clean upholstery, this YA-02E plus ophthalmic operation stretcher is ideally suited for use in ophthalmic clinics and operating rooms.

Main features

- Electric braking system: easy brake locking, orientation and release for convenient cleaning and relocation
- Position memory function with one-touch reset to recall saved stretcher positions
- Retractable push handles at head and leg ends; quick-direction buttons synchronize caster orientation for smooth movement
- Drop-down dual-lever foot controller to adjust stretcher positions and free medical staff's hands
- Head section with longitudinal and angular adjustment for enhanced flexibility
- High-capacity lead-acid battery for long service time
- PU foam headrest and PU artificial leather cushion comfortable and easy to clean
- Molded cast aluminum alloy chassis with high rigidity and elegant appearance

Main details

Headrest

Molded headrest that can be adjusted to various angles both electrically and manually.







Hand remote control

Hand remote control with clearly marked icon buttons and five memory function buttons.



Foot control

The foot-operated control enables hands-free operation and is fitted with an emergency stop button.



Battery

The operating stretcher is equipped with two batteries, enabling an exceptionally long runtime.



Handle

The foldable, stowable handle saves space, while the silver button on the handle controls the direction of the casters.



Stretcher platform

Seamless, three-section PU foam mattress surface for easier cleaning.

Function



Trendelenburg/Reverse Trendelenburg



Back part/leg part



Height adjustment



1 Transport



2 Anesthesia



3 Treatment / Surgery 1



4 Treatment / Surgery 2



5 Flat position



Technical parameter

Model	YA-02E Plus
Overall Size(LWH)	1900×650×550 mm (74.8'' × 25.5'' × 21.6'')
Height Adjustment	550-810 mm (21.6'' - 31.8'')
Head section up/down, electric	25°/-25°
Head section up/down, manual	30°/-30°
Back section up	75°
Leg section down	-75°
Trend/reverse trend	5°/-20°
Sliding	130 mm (5'')
Load capacity	320 kg (706lbs)

ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.

YOUTUBE MEDIK

Pls scan the QR code for more info

info@mediksurgical.com | +86 139 6223 1292

Medik reserves the righto make changes without notice in design,specfcations and models The ony waranty Medik makes is the express writen waranty extended on the sale or rental of its products.

©2026 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 01-EN-2026 EN-US